

MAGNETOM ESPREE



Contents:

Sheet No.	Description
A-101	EQUIPMENT PLAN-DETAILS, LEGEND, AND NOTES
A-102	SAFETY/SERVICE CLEARANCE PLAN
A-501	DETAILS AND NOTES
A-502	DETAILS AND NOTES
S-101	STRUCTURAL PLAN-DETAILS AND NOTES
E-101	ELECTRICAL PLAN(S)-LEGEND AND NOTES
E-102	ELECTRICAL PLAN-LEGEND AND NOTES
E-501	ELECTRICAL DIAGRAMS AND NOTES
M-101	MECHANICAL PLAN AND DETAILS
M-501	CRYOGEN NOTES AND DETAILS

Typical Drawing #: 04103

SIEMENS
SIEMENS MEDICAL SOLUTIONS
51 Valley Stream Parkway
Malvern, PA 19355
www.usa.siemens.com/medical

ARCHITECTURAL NOTES

- 1) ALL PRELIMINARY EQUIPMENT LAYOUTS SUBMITTED BY SIEMENS MEDICAL SOLUTIONS, INC. (SMS HEREAFTER) ARE BASED ON THE RECOMMENDED SPACE NECESSARY FOR THE OPERATION AND SERVICEABILITY OF THE EQUIPMENT BEING PROPOSED. SMS WILL NOT SUBMIT AN EQUIPMENT LAYOUT THAT IS NOT IN THE BEST INTEREST OF BOTH THE CUSTOMER AND SMS. ALL EQUIPMENT LAYOUTS ARE BASED EITHER ON AN ACTUAL SITE LOCATION SURVEY OR ARCHITECTURAL DRAWINGS SUPPLIED TO SMS. SMS WILL NOT BE RESPONSIBLE FOR ANY ALTERATIONS THAT ENCR OACH WITHIN DESIGNATED SAFETY AND SERVICE CLEARANCE ZONES AS INDICATED ON DRAWINGS (IE. PIPE CHASES, VENTILATION DUCTS, CASEWORK, AND SOFFITS, ETC.) MADE BY THE CUSTOMER OR REQUIRED BY A CUSTOMER'S ARCHITECTURAL FIRM ONCE PRELIMINARY DRAWINGS HAVE BEEN SUBMITTED AND APPROVED. DO NOT ALTER ANY SPECIFICATIONS AND/OR DIMENSIONS WITHOUT CONTACTING AND RECEIVING WRITTEN CONFIRMATION FROM SMS PROJECT MANAGER.
- 2) SMS IS NOT AN ARCHITECTURAL OR ENGINEERING FIRM. DRAWINGS SUPPLIED BY SMS ARE NOT CONSTRUCTION DRAWINGS. THEREFORE, THESE DRAWINGS ARE TO BE USED ONLY FOR INFORMATION TO COMPLEMENT ACTUAL CONSTRUCTION DRAWINGS AVAILABLE FROM A CUSTOMER APPOINTED ARCHITECTURAL REPRESENTATIVE OR A CUSTOMER'S ENGINEERING DESIGN GROUP. SMS REQUIRES THAT ONCE THE FINAL CONSTRUCTION DRAWINGS HAVE BEEN PREPARED, THEY SHALL BE MADE AVAILABLE TO SMS PROJECT MANAGER TO VERIFY THAT ALL SMS REQUIREMENTS HAVE BEEN ADHERED TO. THE CUSTOMER'S ARCHITECT AND GENERAL CONTRACTOR SHALL BE ULTIMATELY RESPONSIBLE FOR COMPLIANCE WITH ALL APPLICABLE CODES AND PROFESSIONAL DESIGN REQUIREMENTS.
- 3) THE CUSTOMER IS RESPONSIBLE FOR ALL ROOM AND AREA PREPARATION COSTS, PROFESSIONAL FEES, PERMITS, REPORTS, AND INSPECTION FEES.
- 4) EQUIPMENT WARRANTIES, EXPRESSED OR IMPLIED ON THE PART OF SMS SHALL BE CONTINGENT UPON STRICT COMPLIANCE WITH THE ARCHITECTURAL, STRUCTURAL, ELECTRICAL, MECHANICAL AND RECOMMENDATIONS AND REQUIREMENTS CONTAINED IN THESE DRAWINGS, UNLESS SPECIFIED OTHERWISE.
- 5) ALL DIMENSIONS SHOWN ARE TAKEN FROM FINISHED SURFACES UNLESS SPECIFIED OTHERWISE.
- 6) THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST. ACTUAL PROTECTION REQUIREMENTS SHALL BE SPECIFIED BY A REGISTERED RADIATION PHYSICIST AT CUSTOMER'S ENGAGEMENT AND EXPENSE. RESPONSIBILITY FOR ALL INFORMATION AS TO THE ROOM LOCATION, USE, AND NUMBER OF ANTICIPATED EXAMINATIONS TO BE PERFORMED PER TIME PERIOD SHALL BE PROVIDED TO THE PHYSICIST BY THE CUSTOMER. THE CUSTOMER SHALL FURTHER TAKE ALL RESPONSIBILITY IN THE COMMUNICATION AND COORDINATION OF ACTIVITIES OF THE RADIATION PHYSICIST AND THE ARCHITECTURAL REPRESENTATIVE.
- 7) SMS SHALL BE RESPONSIBLE FOR SMS EQUIPMENT INSTALLATION AND CALIBRATION, CONNECTION AND INSTALLATION OF SMS PROVIDED CABLES, AND CONNECTION OF CONTRACTOR PROVIDED WIRES TO SMS EQUIPMENT. IN THE EVENT THAT SPECIFIC RULES OR LICENSE REQUIREMENTS PROHIBIT THIS, THE CUSTOMER SHALL INITIATE THE SERVICES OF APPROVED OTHER CONTRACTORS AND PAY FOR SELECTED, APPROVED PARTIES TO PERFORM THIS WORK WITH JOB SUPERVISION TO BE PROVIDED BY SMS. CALIBRATION WHEN ACCOMPLISHED OUTSIDE OF NORMAL INSTALLATION SEQUENCES DUE TO CONTRACTOR OR TRADE RULE ACTIONS OR REQUIREMENTS SHALL BE SUPPORTED BY, CHARGED TO, AND ACCEPTED BY THE CUSTOMER AS AN ADDITIONAL INSTALLATION EXPENSE.
- 8) THE CUSTOMER SHALL VERIFY WITH SMS PROJECT MANAGER FINAL INSTALLATION DRAWINGS THE LOCATIONS AND TRAVEL OF ALL ANCILLARY EQUIPMENT TO BE CEILING OR WALL MOUNTED (IE. O.R. LIGHTS, MEDICAL GAS COLUMNS, PHYSIOLOGICAL MONITORING INFUSORS, CRT PLATFORMS, SPRINKLER HEADS, SMOKE DETECTORS, ELECTRICAL OUTLETS, HVAC GRILLES, SPEAKERS, AND GENERAL ROOM LIGHTING, ETC.).
- 9) THE GENERAL CONTRACTOR/CUSTOMER SHALL BE RESPONSIBLE FOR ALL FINAL PAINT, TOUCH-UP AND ANY COSMETIC OR TRIM WORK WHICH NEEDS TO BE OR IS REQUIRED TO BE COMPLETED AFTER THE INSTALLATION OF THE SMS EQUIPMENT AND ANY ASSOCIATED SUPPORT APPARATUS.

SITE READINESS GUIDELINES

- THE FOLLOWING GENERAL CONDITIONS ARE NECESSARY TO HAVE THE STATUS OF "READY SITE":
- 1) PROPER POWER AVAILABLE AT SIEMENS EQUIPMENT POWER CABINET LOCATION AND ALL POWER OUTLETS FUNCTIONING.
 - 2) AIR CONDITIONING/HUMIDIFICATION SYSTEMS COMPLETE, TESTED, AND FUNCTIONING PROPERLY ACCORDING TO SIEMENS SPECIFICATIONS.
 - 3) PROPER LIGHTING INSTALLED AND FUNCTIONING.
 - 4) PLUMBING COMPLETE EXCEPT FOR ANY FINAL CONNECTIONS TO SIEMENS EQUIPMENT.
 - 5) ALL CABLE TRAYS/DUCTS/CONDUITS CORRECTLY SIZED, LOCATED, AND INSTALLED ACCORDING TO THE SIEMENS DRAWINGS.
 - 6) ALL REINFORCEMENT PLATES/UNISTRUT INSTALLED AS REQUIRED.
 - 7) ROOM FOR EQUIPMENT INSTALLATION AND IMMEDIATE VICINITY IS DUST-FREE AND IS TO REMAIN SO FOR THE DURATION OF THE INSTALLATION.
 - 8) A SECURE AREA (APPROXIMATELY 10' x 10') IS AVAILABLE AT EQUIPMENT DELIVERY FOR PARTS AND INSTALLATION TOOLS.
 - 9) CUSTOMER SUPPLIED CAMERAS AND PROCESSORS INSTALLED.
 - 10) CUSTOMER APPROVAL FOR SIEMENS REMOTE SERVICES (SRS) CONNECTION, AND CUSTOMER'S I.T. CONTACT INFORMATION AND IP ADDRESSES ESTABLISHED.
 - 11) WALLS TO BE PRIMED AND PAINTED, FLOORS TO BE TILED EXCEPT IN AREAS OF THE EQUIPMENT BASE PLATES.
- IF THESE CONDITIONS ARE NOT MET, THE SIEMENS PROJECT MANAGER AND THE DESIGNATED SIEMENS INSTALLATION SUPERVISOR SHALL RESCHEDULE THE INSTALLATION START DATE. NOTE: ADDITIONAL COST MAY BE INCURRED BY THE CUSTOMER/CONTRACTOR AND DELIVERY DATES MAY NEED TO BE RESCHEDULED, WHEN THE SIEMENS SITE READINESS GUIDELINES ARE NOT MET.

RESOURCE LIST (SMS USE ONLY)

DESIGNATION	PG NUMBER	DATE
MAGNETOM ESPREE	M6-010.891.01.09.02	2.10

ESPREE
03/01/11

SIEMENS

MAGNETOM ESPREE

TYPICAL FINAL DRAWINGS
MRI SUITE

PROJECT #:

04103

SHEET 1 OF 10
DRAWN BY: B. HERRMANN

SHEET:

A-101

DATE: 03/17/11
CHECKED:

EQUIPMENT LEGEND

NO	DESCRIPTION	SMS SYM	WEIGHT (LBS)	BTU/HR TO AIR	DIMENSIONS (INCHES)			REMARKS
					W	D	H	
1	MRC OPERATING CONSOLE AND KEYBOARD	Ⓐ	132	---	45 11/16	35 1/4	28 3/8	ON CONSOLE/COUNTER
2	COLOR MONITOR FOR MRC	Ⓑ	22	239	18 5/16	16 15/16	4 3/4	
3	HOST PC MRC	Ⓒ	49	2,389	11	27	18 1/8	
4	CONTAINER FOR HOST 500	Ⓓ	238	---	19 5/8	31 1/2	28 3/8	
5	ALARM BOX	Ⓔ	3	---	9	4	9	
6	PATIENT MONITOR (OPTION)	Ⓕ	10	---	13	7	12 1/2	
7	PATIENT SUPERVISION CAMERA (OPTION)	Ⓖ	8	---	3	6 5/8	5 3/4	
8	ESPREE MAGNET WITH COVERS AND PATIENT TABLE	Ⓗ	11,244	7,506	90 5/8	149 5/8	90 5/8	
9	RF-FILTER PLATE	Ⓘ	286	853	46 1/2	35 1/8	21 5/8	
10	MAGNET STOP	Ⓚ	1	---	3	5	3	
11	PATIENT TRANSPORT TROLLEY (OPTION)	Ⓛ	291	---	26 1/2	71 1/2	39 1/2	
12	ELECTRONICS CABINET (GPA, ACC & ACS)	Ⓜ	2,755	13,649	63	25 5/8	77 1/2	TOTAL OF GPA, CCA, CCS
13	SEP CABINET	Ⓝ	750	3,412	25 5/8	25 5/8	73 5/8	
14	POWERWARE 9130 UPS WITH EBM (OPTION)	Ⓖ	186	1,257*	16 7/8	12 7/8	16 1/4	*1,765 ON BATTERIES
15	MULTI MODALITY WORKPLACE COMPUTER (OPTION)	Ⓟ	55	---	10	19 3/4	23 5/8	ON/UNDER COUNTER
16	MULTIMODALITY WORKPLACE MONITOR (OPTION)	Ⓡ	57	---	18	18 1/4	18 3/4	ON CUSTOMER'S COUNTER
17	MEDRAD ICBC INJECTOR STAND AND HEAD (OPTION)	Ⓢ	37.5	---	19 1/4	21 1/2	52 1/4	INJECTOR ON STAND
18	MEDRAD ICBC INJECTOR CRU (OPTION)	Ⓣ	7.9	---	12	10 1/2	10 7/8	ON CUSTOMERS COUNTER
19	MEDRAD ICBC INJECTOR POWER SUPPLY (OPTION)	Ⓤ	3.5	---	10	5	3 1/4	OUTSIDE 5mT FIELD

PROTECTING THE MAGNETIC FIELD

THE SIEMENS MAGNETOM UTILIZES A SUPERCONDUCTIVE MAGNET WITH AN EXTREMELY HOMOGENEOUS FIELD WITHIN THE MAGNET TO PROVIDE DISTORTION-FREE IMAGING. THE PRESENCE OF FERROMAGNETIC MATERIAL WITHIN THE VICINITY OF THE MAGNET CAN ADVERSELY AFFECT THE UNIFORMITY OF THE USEFUL MAGNETIC FIELD. THIS APPLIES TO STATIONARY FERROUS MATERIAL (STRUCTURAL STEEL) WHICH IS TO BE MINIMIZED. STATIONARY STEEL COMPENSATION MAY BE ACHIEVED BY MAGNET POSITIONING AND SELECTIVE USE OF SHIMS. FIELD DISTORTION ENCOUNTERED BY MOVING FERROMAGNETIC OBJECTS IS MORE DIFFICULT TO COMPENSATE AND MAY REQUIRE THE USE OF MAGNETIC SHIELDING.

PROTECTING THE ENVIRONMENT

PROTECTING THE IMMEDIATE ENVIRONMENT FROM THE EFFECT OF THE MAGNETIC FIELD REQUIRES CONSIDERATION. INFORMATION STORED ON MAGNETIC DATA CARRIERS SUCH AS DISKS, TAPES, AND CREDIT CARDS MAY BE ERASED IF IN CLOSE PROXIMITY. CAUTION WITH REGARD TO HEART PACEMAKERS MUST BE EXERCISED. MOST PACEMAKER UNITS EMPLOY A REED RELAY WHICH MAY CHANGE OPERATING MODE WHEN EXPOSED TO AN EXTERNAL MAGNETIC FIELD. THEREFORE, PACEMAKER USERS MUST BE KEPT AT A SPECIFIED DISTANCE FROM THE MAGNET WHICH IS DETERMINED BY THE MAGNETIC FIELD STRENGTH.

MAGNET SITING REQUIREMENTS

IT MUST BE ENSURED THAT THE MAGNET IS LOCATED SO THAT THE STABILITY AND HOMOGENEITY OF THE MAGNETIC FIELD ARE NOT ADVERSELY AFFECTED BY EXTRANEOUS FIELDS AND STATIC OR DYNAMIC FERROMAGNETIC OBJECTS.

X/Y AND Z AXIS	SOURCE OF INTERFERENCE
3'-6"	FLOOR STEEL REINFORCEMENT<20 LBS./ FT ² IRON BEAMS < 66 LBS./FT.
16'-1" / 19'-1"	STRETCHERS UP TO 110 LBS.
13'-2"	A/C CHILLERS
17'-5" / 21'-4"	TRANSPORT DEVICES UP TO 440 LBS.
18'-1 / 24'-8"	VEHICLES UP TO 2,000 LBS.
20'-5" / 29'-7"	ELEVATORS, TRUCKS UP TO 10,000 LBS.
39'-5"/26'-3"	AC TRANSFORMERS LESS THAN 100 KVA
41'-1"/32'-10"	AC TRANSFORMERS LESS THAN 250 KVA
42'-8"/39'-5"	AC TRANSFORMERS LESS THAN 650 KVA
46'-0"/49'-3"	AC TRANSFORMERS LESS THAN 1600 KVA
9'-11"/6'-7"	AC CABLES, MOTORS LESS THAN 100 AMPS
23'-0"/9'-11"	AC CABLES, MOTORS LESS THAN 250 AMPS
131'-2"	ELECTRIC RAILWAY SYSTEMS

FOR IRON OBJECTS LOCATED UP TO 45' FROM THE Z AXIS, THE DISTANCES FOR THE Z AXIS MUST BE USED.
REDUCTION IS POSSIBLE WITH STEEL SHIELDING.

ENVIRONMENTAL/POWER AUDIT

AS AN INDICATION OF OUR COMMITMENT TO QUALITY, SIEMENS MAY, AT NO COST TO YOUR FACILITY, CHECK THE OPERATING ENVIRONMENT AFTER SYSTEM TURNOVER TO DETERMINE IF THE REQUIREMENTS FOR TEMPERATURE, HUMIDITY, POWER, AND GROUNDING ARE MET AS PER SIEMENS' PUBLISHED SPECIFICATIONS. SIEMENS WILL GENERATE A WRITTEN REPORT DETAILING THE ENVIRONMENTAL AND ELECTRICAL CONDITION OF THE SITE AFTER TURNOVER AND WILL SHARE THE REPORT WITH YOU. IN THE EVENT WE IDENTIFY ANY ENVIRONMENTAL/POWER DEFICIENCIES AT THE SITE, YOUR FACILITY WILL BE REQUESTED TO CORRECT DEFICIENCIES WITHIN THIRTY (30) DAYS. SHOULD ANY CORRECTIVE ACTIONS BE NECESSARY, AND UPON REQUEST, SIEMENS WILL PROVIDE GUIDANCE IN AN EFFORT TO FACILITATE RESOLUTION. PLEASE BE ADVISED THAT AFTER 30 DAYS NOTICE ANY REPAIR OR MAINTENANCE SERVICES NECESSITATED BY SEVERE DEFICIENCIES WILL FALL OUTSIDE YOUR WARRANTY COVERAGE.

CEILING HEIGHTS

MAGNET ROOM 7'-11" TECHNICAL MINIMUM
MAGNET ROOM 8'-2" RECOMMENDED MINIMUM
CONTROL ROOM 6'-11" MINIMUM
EQUIPMENT ROOM 7'-3" MINIMUM

ARCHITECTURAL EQUIPMENT PLAN

SCALE: 1/4" = 1'-0"

CONSTRUCTION REQUIREMENTS

THE CUSTOMER/CONTRACTOR IS RESPONSIBLE FOR SUPPLYING AND INSTALLING ALL CONSTRUCTION MATERIALS INCLUDING ELECTRICAL AND MECHANICAL DEVICES REQUIRED BY SIEMENS SPECIFICATIONS AND TO ENSURE THAT THE MATERIAL USED INSIDE THE RF-SHIELDING IS AS FREE OF FERROMAGNETIC PROPERTIES AS POSSIBLE.

ROOM MEASUREMENTS

ALL ROOM MEASUREMENTS AND ROOM DETAIL SPECIFICATIONS MUST BE VERIFIED ON SITE PRIOR TO BEGINNING ANY CONSTRUCTION WORK.

STATE AGENCY REVIEW

PRIOR TO SIEMENS EQUIPMENT INSTALLATION, APPROVAL OF CONSTRUCTION OR STRUCTURAL MODIFICATIONS UTILIZING X-RAY FOR DIAGNOSTIC OR THERAPEUTIC PURPOSES, MUST BE OBTAINED BY THE CUSTOMER FROM THE APPROPRIATE STATE AGENCY, IF APPLICABLE.

STATIC DISSIPATIVE FLOORING

SIEMENS RECOMMENDS "STATIC DISSIPATIVE" FLOOR COVERING WITH AN ELECTRICAL RESISTANCE OF $\leq 10^9$ OHMS IN ALL AREAS WHERE SIEMENS EQUIPMENT IS INSTALLED. A STATIC DISSIPATIVE FLOOR REDUCES THE RISK OF STATIC ELECTRIC DISCHARGES THAT MAY DAMAGE SENSITIVE EQUIPMENT AND COMPONENTS.

DC LIGHTING

THE MAGNETIC FIELD ADVERSELY AFFECTS THE OPERATING LIFE OF LIGHT BULBS LOCATED IN THE IMMEDIATE VICINITY OF THE MAGNET. THE FLAMENT IN THE BULB OSCILLATES WITH THE FREQUENCY OF THE POWER SUPPLY. DURING SCANNING, IT IS RECOMMENDED THAT LIGHT FIXTURES IN THE VICINITY OF THE MAGNET BE CONNECTED TO A DC VOLTAGE SUPPLY.

CASEWORK & ACCESSORY NOTES

1. ALL CASEWORK IS EITHER EXISTING OR IS TO BE DESIGNED, DETAILED, FURNISHED AND INSTALLED BY THE CUSTOMER AND/OR CONTRACTOR. FOLLOW DESIGN RECOMMENDATIONS INCLUDED HEREWITH, AS THEY ARE ESSENTIAL FOR THE SUCCESSFUL INSTALLATION & OPERATION OF THE SIEMENS EQUIPMENT.
2. ALL FURNITURE (CHAIRS, ETC.) FOR THE CONTROL ROOM ARE TO BE PROVIDED BY THE CUSTOMER.

TRANSPORTING REQUIREMENTS

LARGEST ITEM WITHOUT PACKING MATERIAL: MAGNET-11,244 POUNDS
MAGNET AS DELIVERED FROM FACTORY WITHOUT TRANSPORT DEVICE: 7'-4" H. (WITHOUT 90° ELBOW MOUNTED) x 7'-7" W. x 8'-10" L.

STANDARD ROOF OPENING - 9'-2" x 7'-11"

IF TRANSPORTING THE MAGNET UP A RAMP, A 15' MAXIMUM ANGLE MUST BE MAINTAINED.

TO TRANSPORT THE GPA/ACC CABINET (63" x 27" x 78" HIGH; 3307 POUNDS), A MINIMUM ROOM HEIGHT OF 6'-9" WITH TRANSPORT ROLLERS, OR 6'-5" WITHOUT ROLLERS IS REQUIRED.

NOISE LEVELS

SYSTEM ROOM	NOISE LEVEL / dB(A)
CONTROL ROOM	</= 55 (AVERAGE VALUE)
EXAMINATION ROOM	</= 85.4 (8 HOUR AVERAGE) (+3dB(A) TOLERANCE = 92.9dB(A)
EQUIPMENT ROOM	</= 65 (AVERAGE VALUE)

THE PHYSICAL CHARACTERISTICS OF THE MR SYSTEM GENERATE A CERTAIN AMOUNT OF NOISE. THIS TABLE HAS INFORMATION TO INSTALL NOISE ATTENUATION TO MEET ANY STATE/LOCAL/OSHA CODES.

ATTENTION:

—THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.
—THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

—IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.


—ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.
—THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

—ISSUE BLOCK—

SYM DATE DESCRIPTION


SCALE: AS NOTED

REF. #: ---




SAFETY/SERVICE CLEARANCE PLAN

SCALE: 1/4" = 1'-0"




MAGNETOM ESPREE TYPICAL SIDE ELEVATION

SCALE: 1/4"=1'-0"



MAGNETOM ESPREE MAGNET ELEVATIONS

SCALE: NONE



MAGNETOM ESPREE TYPICAL FRONT ELEVATION

SCALE: 1/4"=1'-0"

CEILING HEIGHTS

MAGNET ROOM 7'-11" TECHNICAL MINIMUM
MAGNET ROOM 8'-2" RECOMMENDED MINIMUM
CONTROL ROOM 6'-11" MINIMUM
EQUIPMENT ROOM 7'-3" MINIMUM

ATTENTION:

—THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.
—THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

—IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

—ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.
—THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

ESPREE
03/01/11

SIEMENS
MAGNETOM ESPREE
TYPICAL FINAL DRAWINGS
MRI SUITE



THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.
ALL RIGHTS ARE RESERVED.

PROJECT #:
04103
SHEET 2 OF 10
DATE: 03/17/11

DRAWN BY:
B. HERRMANN
CHECKED:

SHEET:
A-102

SYMBOL	DATE	DESCRIPTION
△		
—ISSUE BLOCK—		
SCALE: AS NOTED	REF. #: ---	

RF DOOR OPENING

IN THE EVENT OF A CATASTROPHIC FAILURE OF THE QUENCH VENT DURING A QUENCH, PRESSURE BUILT UP MAY PREVENT OPENING A DOOR THAT OPENS INTO THE RF ROOM, PREVENTING EVACUATION FROM LIFE THREATENING CONDITIONS.

FOR THIS REASON THE RF DOOR SHOULD OPEN TO THE OUTSIDE OF THE RF ROOM. IF THE DOOR CANNOT OPEN OUT FROM THE RF ROOM, OTHER APPROPRIATE MEANS HAVE TO BE PROVIDED SO THAT THE RF ROOM DOOR IS NOT PREVENTED FROM OPENING DUE TO PRESSURE.

IF THE DOOR OPENS INTO THE RF ROOM, A 24"x24" OPENING FOR PRESSURE EQUALIZATION INTO THE RF ROOM MUST BE INSTALLED. THIS IS MANDATORY. THIS IS NOT AN ESCAPE HATCH. THE PURPOSE OF THE OPENING IS TO RELIEVE PRESSURE AND ALLOW THE MAIN DOOR TO BE OPENED SO THAT OCCUPANTS CAN BE EVACUATED.

THE OPENINGS WILL HAVE PANELS INSTALLED IN THE RF ROOM OR THE DOOR THAT CAN BE UNLOCKED AND OPENED TO THE OUTSIDE IN CASE OF EMERGENCY. THESE PANELS REQUIRE AN RF SEALED

MEASURE AT LEAST 24"x24". WHEN USING RECTANGULAR PANELS, THE SHORTER SIDE SHOULD MEASURE OF MINIMUM OF 24".


TO ENSURE UNOBSTRUCTED VENTING, THIS OPENING CANNOT BE SUBDIVIDED. THIS MEANS THAT, FOR EXAMPLE, RF SEALED HONEYCOMB GRIDS ARE NOT PERMITTED.

EASY REMOVAL OF THE PANEL BY A PERSON HAS TO BE ENSURED AND A MINIMUM DISTANCE OF 40" TO A FIXED OBJECT MUST BE MAINTAINED. THE PANEL SHOULD BE INSTALLED IN AN ACCESSIBLE LOCATION AND ALLOW ESCAPE OF THE LOW DENSITY HELIUM.

AS AN ALTERNATIVE TO AN OUTSWING DOOR, THE STATIONARY OBSERVATION WINDOW IS REPLACED BY A WINDOW OPENING INTO THE CONTROL AREA OR THE DOOR IS REPLACED WITH AN RF SEALED SLIDING DOOR. IT SHOULD BE ENSURED THAT THE DOOR CLOSES IN A WAY THAT ALLOWS IT TO MOVE AWAY FROM THE FRAME IN CASE OF OVERPRESSURE.

F THE DOOR OPENS TO THE OUTSIDE, THE OPENING IN THE RF ROOM IS STILL RECOMMENDED.

THE RF ROOM MANUFACTURER CAN PROVIDE YOU WITH ADDITIONAL RF SEALED ROOM OPENINGS THAT LEAD DIRECTLY TO THE OUTSIDE. HOWEVER, THESE OPENINGS ARE ALSO CONDUITS FOR NOISE GENERATED OUTSIDE THE RF ROOM. UNOBSTRUCTED FLOW THROUGH THIS PIPE MUST BE GUARANTEED.










SAFETY ASPECTS FOR THE RF ROOM:


IT MUST BE POSSIBLE TO LOCK THE RF ROOM (EXAMINATION ROOM) DOOR FROM THE OUTSIDE. IT MUST ALSO BE POSSIBLE TO OPEN THE DOOR FROM THE INSIDE WITHOUT A KEY OR ADDITIONAL DEVICE.

THE RF DOOR IS AN IMPORTANT COMPONENT FOR GOOD IMAGE QUALITY AS WELL AS SAFETY, THE OWNER/OPERATOR OF THE MR SYSTEM MUST MAINTAIN THE RF ROOM AS INSTRUCTED BY THE RF ROOM MANUFACTURER IN ORDER TO GUARANTEE CORRECT FUNCTION OF THE RF DOOR.

NO FERROMAGNETIC ITEMS CAN BE BROUGHT INTO THE RF ROOM AFTER THE MAGNET HAS BEEN RAMPED UP TO FIELD. MAGNETIC ITEMS WILL BECOME ATTRACTED TO THE MAGNET WITH NO WARNING AND DUE TO THE HIGH MAGNETIC FIELD, WILL BECOME MISSILES.

NOTE: FOR DOORS MOVED BY AN AUXILIARY DRIVES (ELECTRICAL OR PNEUMATIC), MANUAL OPERATION HAS TO BE ENSURED. AN OUTSIDE WINDOW SHOULD BE IN THE VICINITY TO ALLOW VENTING EXHAUSTED GAS TO THE OUTSIDE. THE INTEGRITY OF THE RF SHIELD MUST BE TESTED AFTER REMODELLING.





ELECTRICAL DIMENSION PLAN

SCALE: $1/4" = 1'-0"$

CEILING HEIGHTS	
MAGNET ROOM	7'-11" TECHNICAL MINIMUM
MAGNET ROOM	8'-2" RECOMMENDED MINIMUM
CONTROL ROOM	6'-11 MINIMUM
EQUIPMENT ROOM	7'-3" MINIMUM

			<div>SIEMENS</div> <div>MAGNETOM ESPREE</div> <div>TYPICAL FINAL DRAWINGS</div> <div>MRI SUITE</div>					
<div>△</div>			THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.		PROJECT #:	SHEET:		
			ALL RIGHTS ARE RESERVED.		04103			
	SYM	DATE	DESCRIPTION	SCALE:	REF. #:		SHEET 7 OF 10	DRAWN BY:
—ISSUE BLOCK—			AS NOTED	---	DATE:	03/17/11	B. HERRMANN	E-102
					CHECKED:			

ATTENTION:

THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.

THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

ESPREE
03/01/11

Z
O
T
F
S
Z
C
r
O
E
T
Z
-
C
O
D
W
E
r
r

GENERAL NOTES


- IN ORDER TO MAINTAIN THE INTEGRITY OF THE RF-SHIELDED EXAM ROOM, THE GROUND MUST ORIGINATE AT THE SIEMENS ELECTRONICS CABINET GROUND BUS, AND SHOULD BE CONNECTED AT BOTH ENDS UNDER SIEMENS SUPERVISION.
- THE CUSTOMER/CONTRACTOR IS TO FURNISH A FLEXIBLE DEDICATED AND ISOLATED GROUND (GRD) CONDUCTOR FROM THE MAIN FACILITY DISTRIBUTION SOURCE TO THE SIEMENS ELECTRONICS CABINET (ACC). THE GROUND CONDUCTOR MUST BE SIZE AND TYPE TO MATCH THE PHASE CONDUCTORS.
- THE INTERNAL LINE IMPEDANCE IS 100mΩHMS MEASURED AT THE SIEMENS ELECTRONICS CABINET (ACC). ALL WIRING AND CONDUITS ARE TO BE SIZED ACCORDINGLY TO MEET THIS SPECIFICATION. ALL CONDUCTORS ARE TO BE THE SAME SIZE.
- THE EPO (EMERGENCY POWER OFF) MUST PROVIDE REMOTE "EMERGENCY OFF" CONTROL OF SYSTEM POWER. THE EPO'S MUST BE INSTALLED BY A QUALIFIED ELECTRICAL CONTRACTOR ACCORDING TO NATIONAL ELECTRICAL CODE, STATE AND LOCAL REGULATIONS. MEASURES SHOULD BE TAKEN TO DESIGN THE CIRCUIT IN SUCH A WAY THAT IT WILL ALWAYS WORK WHEN THE MEDICAL EQUIPMENT IS POWERED. IT IS RECOMMENDED THAT EPO CIRCUIT BE A NORMALLY CLOSED ENERGIZED CIRCUIT DESIGN AND THAT IT DERIVES ITS SOURCE FROM THE EQUIPMENT INCOMING FEEDER OR IS A DEDICATED BRANCH CIRCUIT, PREFERABLY FROM AN EVERGENCY POWER SOURCE WITH CIRCUIT BREAKER LABELED AND "LOCKED-ON". THE EPO CONFIGURATION DEPICTED IN THIS DRAWING IS ONE EXAMPLE OF A POSSIBLE EPO CONFIGURATION THAT SATISFIES THESE REQUIREMENTS. HOWEVER, THE FACILITY/CONTRACTOR IS SOLELY RESPONSIBLE FOR THE IMPLEMENTATION OF THE EPO AND MUST MAKE THE FINAL DETERMINATION CONSIDERING ALL SITE CONDITIONS AND REGULATORY FACTORS.

THE SECOND SET OF NORMALLY CLOSED CONTACTS ARE TO BE DRY TYPE AND CONNECTED IN SERIES FOR THE UPS. (THEORY OF CIRCUIT) WHEN ANY EPO IS PUSHED, THE NORMALLY CLOSED UPS CONTACT WILL OPEN SENDING A SIGNAL TO THE UPS SHUTTING THE UPS DOWN.

LEGEND

- A E1SED43B125 - BREAKER AND ENCLOSURE
B U01FD60 - 120VAC UNDERVOLTAGE TRIP
D KT8050P - 480 VOLTS TO 120 VOLTS STEP DOWN TRANSFORMER FOR UNDERVOLTAGE TRIP
G GROUND BAR
- EMT CONDUIT
- - - - - ALUMINUM CONDUIT
- · - · - · FLEXIBLE CONDUIT WITH DIELECTRIC AT ACC

RECOMMENDED PART NUMBERS SHOWN ARE SIEMENS ENERGY AND AUTOMATION (PURCHASED FROM YOU LOCAL DISTRIBUTOR) OR EQUIVALENT. TO LOCATE A DISTRIBUTOR, VISIT THIS WEBSITE-[HTTP://WWW.SEA.SIEMENS.COM](http://www.sea.siemens.com). ALL PARTS TO BE PROVIDED BY CUSTOMER/CONTRACTOR.



POWER AND GROUNDING REQUIREMENTS

SCALE: NONE

SIEMENS REMOTE SERVICES (SRS)

TO ENSURE THE UPTIME OF YOUR SYSTEM DURING THE WARRANTY PERIOD (AND BEYOND WITH A SERVICE AGREEMENT), SIEMENS REMOTE SERVICES (SRS) REQUIRES REMOTE LOCAL AREA NETWORK ACCESS TO SIEMENS SYSTEMS.

SRS REQUIRES ONE OF THE FOLLOWING CONNECTION METHODS:


(PREFERRED) VPN CONNECTION

THE PREFERRED CONNECTION METHOD IS (VPN) VIRTUAL PRIVATE NETWORK (WHERE THE CUSTOMER HAS AVAILABLE A VPN CAPABLE FIREWALL OR OTHER VPN APPLIANCE). THIS METHOD PROVIDES THE POSSIBILITY FOR REMOTE SYSTEM DIAGNOSTICS WITHOUT ADDITIONAL HARDWARE. PLEASE CONTACT SIEMENS REMOTE SERVICES (800-888-SIEM) TO DETERMINE IF THIS METHOD IS SUITABLE FOR YOUR SITE.

(OPTIONAL) SRS ROUTER CONNECTION

- THE SRS ROUTER IS SUPPLIED BY SIEMENS AND INSTALLED AT THE CUSTOMER'S SITE, WHILE STILL REMAINING THE PROPERTY OF SIEMENS. THE CUSTOMER'S NETWORK ADMINISTRATOR AND SIEMENS REMOTE SERVICES SHALL DETERMINE THE TYPE AND LOCATION OF THE SRS ROUTER REQUIRED.
- THE SRS ROUTER IS CONNECTED TO AN ANALOG MODEM THAT IS SUPPLIED BY SIEMENS, WHICH THEN IN TURN IS CONNECTED TO AN ANALOG PHONE LINE THAT IS SUPPLIED BY THE CUSTOMER. ONE SRS ROUTER ALLOWS REMOTE DIAGNOSTICS TO MULTIPLE MEDICAL SYSTEMS.
- THE SRS ROUTER SHOULD BE INSTALLED IN A SECURE LOCATION (CUSTOMER'S NETWORK COMPUTER ROOM) THAT HAS LIMITED ACCESS. IT CAN BE LOCATED ON A SHELF, TABLE, OR IN A CABINET. THE CONNECTION CABLES (WITH INDICATED LENGTHS BELOW) ARE INCLUDED WITH DELIVERY.

SRS ROUTER CONNECTION DIAGRAM



NOTE: ALL POWER OUTLETS ARE SUPPLIED/INSTALLED BY CUSTOMER.

- ETHERNET SWITCH OR HUB, SUPPLIED BY CUSTOMER
- SRS ROUTER, SUPPLIED BY SIEMENS (SIZE: 11.2"W x 8.7"D x 5.5"H, WEIGHT: 2 LBS.)
- ANALOG MODEM, SUPPLIED BY SIEMENS
- ANALOG PHONE LINE, SUPPLIED BY CUSTOMER

SIEMENS REMOTE SERVICE

SCALE: NONE

LIGHTING GUIDELINES

EXAM, CONTROL, AND EVALUATION ROOMS:

- THE ROOM LIGHTING MUST REMAIN FUNCTIONAL WHEN THE MR SYSTEM IS SWITCHED OFF AND/OR WHEN EMERGENCY SHUTDOWN BUTTONS ARE ACTIVATED.
- IT MUST BE POSSIBLE TO CONTROL THE INTENSITY OF ILLUMINATION OF APPROXIMATELY 46 FOOT-CANDLES THROUGH GROUP CONNECTION. ALL LIGHTS IN THE EXAMINATION ROOM ARE CONNECTED TO A COMMON SWITCH IN THE CONTROL ROOM. THERE SHOULD BE SEPARATE SWITCHES IN THE EXAMINATION ROOM FOR THE GROUPS OF LIGHTS ABOVE AND NEAR THE PATIENT TABLE, AS WELL AS FOR THE GROUP OF LIGHTS ABOVE THE MAGNET. THE LOCATION OF THESE SWITCHES IS AT THE OWNER'S DISCRETION.

EXAMINATION ROOM:

- THE MAGNETIC FIELD ADVERSELY AFFECTS THE OPERATING LIFE OF LIGHT BULBS LOCATED IN THE IMMEDIATE VICINITY OF THE MAGNET. THE FILAMENT IN THE BULB OSCILLATES WITH THE FREQUENCY OF THE POWER SUPPLY. DURING SCANNING, IT IS RECOMMENDED THAT LIGHT FIXTURES IN THE VICINITY OF THE MAGNET (IN EXAMINATION ROOM) BE CONNECTED TO A DC VOLTAGE SUPPLY. THE RESIDUAL RIPPLE OF THE DIRECT VOLTAGE SHOULD BE ±5%, WHEN INSTALLING THE LIGHT SOCKET, ENSURE THAT THE POLARITY IS CORRECT.
- FLUORESCENT LIGHTS, ENERGY- SAVING LIGHTS, AND DIMMERS ARE NOT PERMITTED.

CONTROL AND EVALUATION ROOM:

FOLLOW THE APPROPRIATE GUIDELINES FOR LIGHTING IN ROOMS WITH MONITOR WORKSTATIONS.

EQUIPMENT ROOM:

ILLUMINATION INTENSITY SHOULD BE APPROXIMATELY 46 FOOT-CANDLES.

ELECTRICAL INSTALLATION NOTES

- INSTALL THE MR SYSTEM CIRCUIT BREAKER NEAR OR IN THE EQUIPMENT ROOM. THE PERMITTED FRINGE FIELD FOR THE PANEL IS UP TO 3mT. IF THE FRINGE FIELDS HAVE HIGHER VALUES, MAGNETIC SHIELDING MUST BE PROVIDED OR THE DISTANCE FROM THE MAGNET MUST BE INCREASED.
- A UNDERVOLTAGE TRIP CIRCUIT BREAKER FOR SWITCHING MAIN POWER AND OFF SHOULD BE INSTALLED IN THE MAIN BREAKER PANEL. INSTALL EMERGENCY SHUTDOWN BUTTONS IN INDIVIDUAL ROOMS ACCORDING TO THE ELECTRICAL INSTALLATION PLAN. THE BUTTONS INTEGRATED IN THE ALARM BOX (AB) IN THE CONTROL ROOM SHOULD BE USED FOR SWITCHING THE SYSTEM ON AND OFF.
- WALL RECEPTACLES MADE OF FERROMAGNETIC MATERIALS ARE NOT PERMITTED IN THE EXAM ROOM. PERIPHERAL UNITS (SUCH AS SERVO-VENTILATORS) NOT APPROVED FOR USE IN A HIGH MAGNETIC FIELD ENVIRONMENT CAN INFLUENCE THE MAGNETIC FIELD, COMPROMISING IMAGE QUALITY. THE CUSTOMER IS RESPONSIBLE FOR INSTALLATION AND USE OF RECEPTACLES IN THE EXAM ROOM. INSTALLATION OF RECEPTACLES AND THE FILTERS REQUIRED ARE TO BE COORDINATED WITH THE RF SHIELDING SUPPLIER.
- THE EMERGENCY POWER OFF BUTTONS (EPO) ARE TO BE SIEMENS 52PA2W2A OR EQUIVALENT "EMERGENCY STOP" MUSHROOM BUTTONS WITH PUSH LOCK AND PULL TO RELEASE.
- THE RF SHIELD MUST BE FITTED WITH A GROUND STUD OR BUS BAR, LOCATED WITHIN 24" OF THE AUXILIARY FILTERS FOR ROOM LIGHTS AND OUTLETS, SUPPLIED AND INSTALLED BY THE RF SHIELD SUPPLIER.
- ALL CUSTOMER OR CUSTOMER/CONTRACTOR SUPPLIED AC POWER ENTERING THE EXAMINATION ROOM (I.E. OUTLETS, EPO, ETC.) SHOULD BE SUPPLIED VIA AN ISOLATION TRANSFORMER. THE ISOLATION TRANSFORMER SECONDARY WINDING GROUND CONDUCTOR SHOULD BE CONNECTED TO THE RF SHIELD GROUND STUD OR BUS BAR. SEE NOTES 5 ABOVE AND THE AUXILIARY AC POWER FOR EXAMINATION EXAMINATION ROOM DETAIL.

CONDUITS AND RACEWAYS

- ALL POWER CONDUCTORS SUPPLIED BY THE CUSTOMER/CONTRACTOR SHALL BE INSTALLED IN METAL RACEWAY. 600 VOLT CLASS, STRANDED TYPE THIN-WALL, RATED FOR 75°C (165°F) OPERATION. RECOMMEND MINIMUM 5 FEET WIRE TAILS AT ALL OUTLET POINTS WITH WIRE IDENTIFICATION TAGGED AT BOTH ENDS FOR FINAL CONNECTION BY SIEMENS MEDICAL SYSTEMS.
- THE CABLE GROUPS INCLUDED WITH THE MAGNETOM SYSTEM MAY BE ROUTED IN THE SAME CABLE TRAY IF PROVIDED WITH AN 8" SEPARATION BETWEEN SMALL SIGNAL LINES, GRADIENT CABLES, AND THE RF TRANSMIT CABLE. A 24" WIDE LADDER TYPE CABLE TRAY IS RECOMMENDED. CABLES SHOULD NOT BE BUNDLED TOGETHER.
- NOTE THE CABLE CONNECTOR SIZES (LARGEST CONNECTOR SIZE IS 2 1/2" x 2 1/2") FOR CABLE FEED-THROUGHS AND CABLE DUCTS.
- THE CABLE LENGTHS SPECIFIED ARE THE STANDARD LENGTHS.
- THE SIEMENS SYSTEM CABLES ARE NOT PLENUM RATED AND SHOULD NOT BE RUN UNPROTECTED IN AN AIR PLENUM UNLESS ENCLOSED IN A SEALED CABLE TRAY OR CONDUIT.

POWER QUALITY NOTES

- IT IS THE CUSTOMER'S RESPONSIBILITY TO COMPLY WITH THE POWER QUALITY REQUIREMENTS FOR SIEMENS MEDICAL SYSTEMS EQUIPMENT.
- THE ELECTRICAL FEEDER TO THE SIEMENS MEDICAL SYSTEM EQUIPMENT MUST FEED ONLY THE IMAGING SYSTEM AND BE KEPT SEPARATE FROM ELECTRICAL FEEDERS TO HVAC, MOTORS, PUMPS, COMPRESSORS, ELEVATORS, AND OTHER POTENTIAL SOURCES OF ELECTRICAL INTERFERENCE.
- THE ELECTRICAL FEEDER TO THE IMAGING SYSTEM MUST BE RUN DIRECTLY TO A MAIN FACILITY DISTRIBUTION PANEL OR TO THE FACILITY SERVICE ENTRANCE, WITH NO OTHER LOADS POWERED FROM THIS FEEDER.
- IN ORDER TO COMPLY WITH IMAGING SYSTEM POWER QUALITY REQUIREMENTS, ADDITIONAL POWER CONDITIONING DEVICES MAY BE REQUIRED. EXAMPLES INCLUDE VOLTAGE REGULATORS, TRANSFORMERS, SURGE PROTECTIVE DEVICES, FILTERS, AND/OR UNINTERRUPTIBLE POWER SUPPLIES (UPS). RECOMMENDATIONS FOR THE INSTALLATION OF ELECTRONIC EQUIPMENT CAN BE FOUND IN IEEE STANDARD 1100-1999 "POWERING AND GROUNDING ELECTRONIC EQUIPMENT."
- POWER CONDITIONING DEVICES NOT APPROVED BY SIEMENS MEDICAL SYSTEMS MAY NOT BE COMPATIBLE WITH THE MAGNETOM SYSTEM. "FERRORESONANT" POWER CONDITIONING EQUIPMENT RE-APPLIED FROM PREVIOUS GENERATION SYSTEMS IS ALSO GENERALLY EXCLUDED DUE TO HIGHER POWER REQUIREMENTS OF THE NEWER SYSTEMS.

GROUNDING NOTES


EQUIPMENT GROUND CONDUCTOR TO COMPLY WITH THE FOLLOWING:

- SIZED EQUIVALENT TO THE PHASE CONDUCTORS (FULL SIZE GROUND).
- DERIVED FROM THE ELECTRICAL SERVICE, TRANSFORMER OR MAIN DISTRIBUTION PANEL FEEDING THE SIEMENS EQUIPMENT.
- RUN IN THE SAME CONDUIT, TROUGH OR RACEWAY AS THE PHASE CONDUCTORS.
- CONTINUOUS, WITH NO BREAKS OR USE OF CONDUIT, CHASSIS OR EARTH AS THE SOLE GROUNDING PATH.
- BONDED TO CHASSIS AND/OR CONDUIT IN ACCORDANCE WITH THE NEC REQUIREMENTS.
- MINIMIZE CONNECTIONS OR TERMINALS TO ENSURE CONTINUITY OVER THE LIFE OF THE INSTALLATION.
- AS A NORM, THERE SHOULD NOT BE ANY CURRENT PRESENCE ON THE GROUND CONDUCTOR, BUT IT IS ACCEPTABLE TO HAVE <500mA DURING OPERATION OF THE IMAGING EQUIPMENT.
- THERE MAY BE SOME APPLICATIONS WHICH REQUIRE AN ISOLATED GROUND AS PER NEC 250-96B.


MR GROUNDING NOTES

THE INTERNAL GROUND WIRING OF THE MR SYSTEM MUST BE INSTALLED WITH MINIMUM GROUND LOOPS. THIS IS TO PREVENT NOISE CURRENTS AND GENERAL DISTURBANCES FROM FLOWING THROUGH THE GROUNDING PATH.

TO ACHIEVE SUCH GROUNDING, THREE MAJOR GROUND POINTS SHOULD BE USED.



THE PROPER ROUTING OF CABLES IS ESSENTIAL TO ACHIEVE GOOD IMAGE QUALITY. RF CABLES MUST BE SEPARATED FROM FIBER OPTIC BY AT LEAST 12" AND FROM THE GRADIENT CABLES BY AT LEAST 12". FIBER OPTIC CABLES MUST ALSO BE SEPARATED FROM THE GRADIENT CABLES BY AT LEAST 12". THIS SHOWS RACEWAY/CABLE ROUTING.



CABLE DESIGNATIONS ARE SHOWN AS AN EXAMPLE. ANY CATEGORY CABLE CAN BE LOCATED IN ANY OF THE COMPARTMENTS OF THE RACEWAY AS LONG AS CORRECT SEPARATIONS ARE MAINTAINED.

WHEN ROUTING RACEWAYS, DO NOT EXCEED THE MAXIMUM LENGTHS LISTED IN DETAIL E-501/2. EXCESS CABLE SHOULD BE ROUTED IN THE RACEWAY IN A MEANDERING METHOD, NEVER ROLLED IN LOOPS.


THE BENDING RADIUS FOR THE CABLES MUST BE MAINTAINED.
TRANSMITTER CABLE - 5" WHEN BENT ONCE.
TRANSMITTER CABLE - 14.25 WHEN BENT SEVERAL TIMES.
FIBER OPTIC CABLE - 6"
GRADIENT CABLE - 5.5" (ONLY WITH EXTENDED CABLE SET)
FIBER OPTIC CABLE FOR PATIENT OBSERVATION - 2"

CABLE SEPARATION

SCALE: NONE

AUXILIARY AC POWER FOR EXAMINATION ROOM

SCALE: NONE



ATTENTION:

THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.
THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.
THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

ESPREE
03/01/11

SIEMENS

MAGNETOM ESPREE

TYPICAL FINAL DRAWINGS
MRI SUITE

PROJECT #:

04103

SHEET 8 OF 10

DRAWN BY: B. HERRMANN

CHECKED:

SHEET:

E-501

SYMBOL

DATE

DESCRIPTION

ISSUE BLOCK


SCALE: AS NOTED

REF. #:

DATE: 03/17/11

DATE: 03/17/11

CHECKED:



MECHANICAL PLAN

SCALE: 1/4" = 1'-0"

CHILLED WATER REQUIREMENTS

WATER REQUIREMENTS TO BE MEASURED AT THE SEP CABINET.	
FLOW RATE:	23.78-29.05 GPM
WATER TEMPERATURE:	48°F ±4°F
BTU DISCHARGE TO THE WATER	163,793 BTU/HR
WATER PRESSURE	MAXIMUM 87 PSI
LOSS OF PRESSURE FOR SEP CABINET	14.5 PSI MAXIMUM
CHILLED WATER ACIDITY RANGE	6 pH TO 8 pH
CHILLED WATER HARDNESS	<250 ppm CALCIUM CARBONATE
CHLORINE GAS CONCENTRATION	<200 ppm
FILTRATION	500 µm

FOR INSTALLATION OF A KRAUS KCC 215 CHILLER, IT IS THE RESPONSIBILITY OF THE CUSTOMER/MECHANICAL CONTRACTOR TO FLUSH PROVIDE A MIXTURE OF WATER WITH 33%-38% ETHYLENE GLYCOL PRIOR TO CHILLER START UP. DO NOT USE PROPYLENE GLYCOL OR AUTOMOTIVE ANTI-FREEZE.

THE AMOUNT OF THE MIXTURE MUST FILL THE CHILLER, MR SYSTEM AND PIPING (SUPPLY AND RETURN), SEE EXAMPLES BELOW.

(1) GALLON OF UNDILUTED GLYCOL, OR (2) GALLONS OF WATER/GLYCOL MIXTURE MUST REMAIN ON SITE FOR USE AFTER START UP.

MIXTURE VOLUME INCLUDING SUPPLY & RETURN+15 GAL. CHILLER & MR			
PIPE DIAMETER	TOTAL LENGTH	MIXTURE VOLUME	GLYCOL NEEDED
2"	100'	31.3 GALLONS	11.9 GALLONS
2"	200'	47.6 GALLONS	18.1 GALLONS
2.5"	100'	40.5 GALLONS	15.4 GALLONS
2.5"	200'	66.0 GALLONS	25.1 GALLONS

MIXTURE VOLUME = 3.14 x (PIPE RADIUS)² x PIPE LENGTH + 15 GALLONS. GLYCOL AMOUNT = 35-38% OF MIXTURE VOLUME.

CHILLED WATER SUPPLY

A CHILLED WATER SUPPLY IS REQUIRED TO THE MRI SYSTEM 24 HOURS A DAY, YEAR ROUND FOR HTE COLD HEAD AND GRADIENT SYSTEMS. THIS CAN BE PROVIDED BY A CENTRAL CHILLED WATER SUPPLY OR A SEPARATE STAND ALONE CHILLER THAT MEETS THE STATED REQUIREMENTS. THE CHILLED WATER CAN ALSO BE SUPPLIED BY A DEDICATED KRAUS KSC 215 CHILLER AND INTERFACE PANEL.

WITHOUT THE USE OF A DEDICATED KRAUS CHILLER, A SEP (SYSTEM SEPARATOR CABINET), MUST BE INCLUDED WITH THE SIEMENS ORDER. THE PIPE SIZE BETWEEN THE KRAUS CHILLER AND INTERFACE PANEL, OR BETWEEN THE WATER SUPPLY AND SEP MUST BE 2 INCH UP TO 82 FEET, 2-1/2 INCH UP TO 148 FEET, CONSULT FOR LONGER PIPE. PERMISSIBLE MATERIALS THAT CAN BE USED FOR THE PIPING ARE: STAINLESS STEEL (V2A, V4A), NON-FERROUS METAL (COPPER, BRASS), SYNTHETIC MATERIAL, PLASTICS, BRAZING SOLDER, HARD SOLDER, OR FITTING SOLDER TYPE 3 AND 4. THERE ARE MATERIALS THAT MAY CAUSE DAMAGE TO THE COOLING SYSTEM AND CANNOT BE USED, THESE MATERIALS ARE ALUMINUM, IRON, CARBON STEEL, ZINC, ZINC PLATED STEEL, OR STANDARD STEEL PIPES.

THESE REQUIREMENTS ARE REQUIRED FOR NEW INSTALLATIONS, IF EXISTING WATER PIPES COMPLY WITH SIEMENS WATER SPECIFICATIONS, THEY DO NOT NEED TO BE REPLACED.

NORMAL TAP WATER MUST BE AVAILABLE FOR FILLING THE SECONDARY WATER CIRCUIT. THERE SHALL BE A HOSE BIB LOCATED WITHIN 65' OF THE SEP, IFP, ACC OR THE KRAUS CHILLER.

THE SUPPLY AND RETURN CHILLED WATER PIPES MUST BE LABELED. THE LOCATION OF THE LABELS MUST BE AT ALL CONNECTION AND REFILLING POINTS AND MUST CONTAIN FLOW DIRECTION AND CONTENTS.

ENVIRONMENTAL REQUIREMENTS

- AIR CONDITIONING IS TO PROVIDE A TEMPERATURE OF 70°F ±5°F IN THE EXAM ROOM, 70°F±10°F IN THE EQUIPMENT & CONTROL AREAS, RELATIVE HUMIDITY OF 40-60% (NON-CONDENSING) IS REQUIRED EXAMINATION ROOM AND 40-80% (NON-CONDENSING) IN ALL OTHER AREAS WHERE SIEMENS EQUIPMENT IS INSTALLED. THESE CONDITIONS ARE TO BE MET AT ALL TIMES: 24 HOURS A DAY, 7 DAYS A WEEK.
- A DEDICATED AIR CONDITIONING AND HUMIDIFICATION SYSTEM IS RECOMMENDED FOR THE EXAM ROOM. A MINIMUM AIR EXCHANGE RATE OF 6 TIMES PER HOUR FOR THE EXAM ROOM IS REQUIRED. IT IS RECOMMENDED TO INSTALL A FRESH AIR SYSTEM WITH 30%-50% FRESH AIR INTAKE.

AIR SUPPLY AND RETURN ABOVE THE FINISHED CEILING IN THE EXAM ROOM IS RECOMMENDED. EACH ROOM SHOULD HAVE A DEDICATED CONTROL AND SENSOR TO MONITOR AND ADJUST THE AIR.


3) THE HEAT INTO THE EXAM ROOM IS LESS THAN 10,236 BTU/HR. THE HEAT INTO THE EQUIPMENT ROOM IS TYPICALLY 8,530 BTU/HR, MAXIMUM 17,060 BTU/HR. THIS HEAT DISSIPATION IS FROM THE SIEMENS EQUIPMENT ONLY. AUXILIARY SUPPORT EQUIPMENT (ie UPS) AND LIGHTING MUST BE CONSIDERED FOR TOTAL HEAT LOADS.

- IT IS IMPORTANT FOR FRESH AIR INTAKE SYSTEMS TO EXHAUST AIR DIRECTLY OUT OF THE BUILDING. THE EXHAUST AIR MUST NOT BE DEFLECTED INTO ANOTHER ROOM. THE MAGNET ROOM EXHAUST AIR SHOULD BE INSTALLED AT LEAST 6'-6" ABOVE FINISHED FLOOR.
- THE AIR INTAKE OF THE AIR CONDITIONING SYSTEM MUST NOT BE LOCATED IN THE VICINITY OF THE QUENCH VENT EXHAUST.

- IF THE INPUT DRAWS UPON AIR FROM OUTSIDE THE BUILDING, IT IS RECOMMENDED TO INSTALL AN ON-SITE FILTER TO REMOVE DUST PARTICLES GREATER THAN 10 MICRONS.

PIPING SCHEMATIC FOR DEDICATED CHILLED WATER

SCALE: NONE



NOTES:


ALL PIPING AND PLUMBING FIXTURES SHALL BE FURNISHED, INSTALLED, PRESSURE TESTED AND CHANGED BY THE MECHANICAL CONTRACTOR PRIOR TO THE DELIVERY AND INSTALLATION OF THE SIEMENS SUPPLIED AND INSTALLED EQUIPMENT UNLESS SPECIFIED OTHERWISE.

AT THE HIGHEST POINT OF THE WATER SUPPLY PIPE FROM THE KRAUS CHILLER AN AUTOMATIC DEAERATION DEVICE (AIR VENT) WITH BALL VALVE MUST BE INSTALLED BY THE MECHANICAL CONTRACTOR. LEAK TEST ALL PIPING WITH A MIXTURE OF R-22 TRACE GAS AND NITROGEN. DO NOT PERFORM LEAK TEST WITH WATER.


THE MECHANICAL ENGINEER OF RECORD SHALL BE ULTIMATELY RESPONSIBLE FOR THE SITE SPECIFIC DESIGN AND SPECIFICATION OF THE MECHANICAL AND PIPING SYSTEMS AS SHOWN AND SHALL BE IN ACCORDANCE WITH ALL APPLICABLE LOCAL, STATE AND NATIONAL CODES. ALL WORK SHALL BE PERFORMED BY THE MECHANICAL CONTRACTOR AND SHALL BE SUBJECT TO COMPLIANCE WITH ALL APPLICABLE LOCAL, STATE AND NATIONAL CODES.

THE SUPPLY AND RETURN PIPES FROM THE CHILLED WATER SUPPLY TO THE SEP/IFP MUST BE LABELED TO SHOW FLOW DIRECTION AND CONTENT (WATER/GLYCOL).

A TAP WATER SUPPLY MUST BE AVAILBE WITHIN 60' OF THE SEP/IFP AND CHILLER CONNECTION FOR FILLING THE CIRCUIT.



ARRANGEMENT WITH FACILITY CHILLED WATER OR STAND ALONE CHILLER



ARRANGEMENT WITH KKT KRAUS KSC 215 CHILLER

MECHANICAL NOTES

- THE AIR H.V.A.C. SYSTEM MUST OPERATE FOR A MINIMUM OF 48 CONSECUTIVE HOURS PRIOR TO THE DELIVERY OF THE EQUIPMENT.
- THE FILTERS MUST BE CHANGED IMMEDIATELY PRIOR TO THE DELIVERY OF THE EQUIPMENT.
- SMS REQUIRES THE USE OF A DEDICATED H.V.A.C. SYSTEM FOR THE EQUIPMENT ROOM TO BE LOCATED, SIZED AND SPECIFIED BY THE MECHANICAL ENGINEER OF RECORD AND TO BE SUPPLIED AND INSTALLED BY THE MECHANICAL CONTRACTOR.
- SMS RECOMMENDS THAT THE CUSTOMER PROVIDE AND INSTALL AN OXYGEN MONITORING SYSTEM WITH VISUAL AND AUDIBLE ALARMS TO INDICATE WHEN THE OXYGEN CONTAINED IN AMBIENT AIR FALLS BELOW PRE-PROGRAMMED SAFETY LEVELS WITH THE SENSOR TO BE LOCATED IN THE SCAN ROOM IN THE AREA DESIGNATED FOR CRYOGEN FILLING.
- THE SIEMENS ACTIVE SHIELDED MAGNET RECIRCULATES LIQUID HELIUM, ELIMINATING THE NEED FOR A DEDICATED CRYOGEN STORAGE AREA. THE RECIRCULATING SYSTEM SIGNIFICANTLY REDUCES THE HELIUM "BOIL OFF". THE MAGNET WILL REQUIRE OCCASIONAL FILLING. A DELIVERY ROUTE FOR CRYOGEN DEWARs MUST BE ESTABLISHED. A MINIMUM 36" CLEARANCE IS REQUIRED.

FIRE CONTROL NOTES

- SIEMENS HAS NO SPECIFIC REQUIREMENT FOR FIRE PROTECTION. FIRE PROTECTION REQUIREMENTS SHALL BE IN ACCORDANCE WITH LOCAL CODES AND CUSTOMER'S INSURANCE REQUIREMENTS. ALL FIRE PROTECTION SYSTEMS SHALL BE DEFINED BY THE ARCHITECT OF RECORD WITH DESIGN, SPECIFICATION AND DETAILING OF THE FIRE PROTECTION SYSTEM BY THE MECHANICAL ENGINEER OF RECORD IN ACCORDANCE WITH SIEMENS GUIDELINES AS STATED HEREIN.

THE ELECTRONIC EQUIPMENT OF THE MR SYSTEMS WILL BE DAMAGED BY WATER. REDUCTION OR ELIMINATION OF WATER USED FOR FIRE SUPPRESSION WILL REDUCE POTENTIAL WATER DAMAGE. PRE-ACTION INERT GAS, OR HALOCARBONS OR OTHER METHODS CAN REDUCE OR ELIMINATE WATER. REFER TO YOUR FIRE PROTECTION PROFESSIONAL.

- THE USE OF SMOKE DETECTORS INSIDE OF THE MR EXAMINATION ROOM IS NOT RECOMMENDED. SMOKE DETECTORS, BY DESIGN, CAN GENERATE NOISE THAT MAY INTERFERE WITH THE MRI EXAMINATION AND CAUSE IMAGE ARTIFACTS. IF THE USE OF A SMOKE DETECTOR IN THE EXAMINATION ROOM IS MANDATED BY LOCAL REQUIREMENTS, SPECIAL NOISE TESTS MUST BE PERFORMED BY SIEMENS SERVICE AFTER THE MRI IS OPERATIONAL. MRI EQUIPMENT PERFORMANCE PROBLEMS DUE TO SMOKE DETECTORS ARE THE RESPONSIBILITY OF THE CUSTOMER AND ARE NOT COVERED UNDER WARRANTY OR SERVICE AGREEMENT.

- ALL MATERIAL USED INSIDE THE MAGNET ROOM SHALL BE NON-MAGNETIC.

- ALL PENETRATIONS IN THE RF CABIN/SHIELD SHALL BE THROUGH A WAVEGUIDE TO BE EQUIPPED WITH A SIEMENS APPROVED DIELECTRIC COUPLER ON BOTH ENDS OF THE WAVEGUIDE. ALL WAVEGUIDES SHALL BE DESIGNED, DETAILED AND SPECIFIED BY THE RF CABIN/SHIELD CONTRACTOR WITH ALL LOCATIONS TO BE DETERMINED BY THE ARCHITECT AND MECHANICAL ENGINEER OF RECORD TO BE ESTABLISHED IN A PRE-PLANNING MEETING PRIOR TO THE DESIGN, SPECIFICATION, AND FABRICATION OF THE RF CABIN/SHIELD.

- EACH ELECTRICAL PENETRATION OF THE RF CABIN/SHIELD FOR ELECTRICAL SERVICING OF THE FIRE PROTECTION SYSTEM SHALL BE THROUGH AN RFI FILTER TO BE SUPPLIED BY THE RF SHIELD CONTRACTOR WITH FILTER LOCATIONS TO BE DETERMINED BY THE ARCHITECT AND THE ELECTRICAL ENGINEER OF RECORD TO BE ESTABLISHED IN A PRE-PLANNING MEETING PRIOR TO THE DESIGN, SPECIFICATION AND FABRICATION OF THE RF CABIN/SHIELD.

- IT IS PERMISSIBLE TO RUN "BLACK PIPE" UP TO THE DIELECTRIC COUPLER ON THE OUTSIDE OF THE RF SHIELD.

- THERE MUST BE NO GROUND CONNECTIONS MADE DURING THE INSTALLATION OF EITHER THE PIPING OR ELECTRICAL FOR THE FIRE PROTECTION SYSTEM.

- THE USE OF HALON IS NOT ACCEPTABLE.

- THE LOCATION OF FIRE CONTROL SYSTEM COMPONENTS SHALL BE COORDINATED THROUGH THE ARCHITECT OF RECORD WITH ALL LOCATIONS TO BE COORDINATED WITH SIEMENS EQUIPMENT LOCATIONS AS SHOWN ON THE 1/4" SCALE EQUIPMENT LOCATION PLAN.

- THE FIRE CONTROL CONTRACTOR SHALL VERIFY EQUIPMENT MOUNTING PROCEDURES AND LOCATIONS ON ANY WALLS CONTAINING RF SHIELDING WITH THE SIEMENS PROJECT MANAGER PRIOR TO THE COMMENCEMENT OF WORK.

ESPREE
03/01/11

SIEMENS

MAGNETOM ESPREE

TYPICAL FINAL DRAWINGS
MRI SUITE

PROJECT #:

04103

SHEET:

M-101

THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.

ALL RIGHTS ARE RESERVED.

SHEET 9 OF 10

DRAWN BY: B. HERRMANN

DATE: 03/17/11

CHECKED:

SCALE: AS NOTED

REF. #: ---

ATTENTION:

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.
- THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.
- THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.


CRYOGEN NOTES

- 1) "CRYOGENS" IS A TERM USED TO IDENTIFY THE REFRIGERANT USED TO MAKE THE MAGNET "SUPER-CONDUCTING". IN THIS APPLICATION, LIQUID AND GASEOUS HELIUM. IN THEIR LIQUID STATE THE CRYOGENS CAN REACH 4.2°K (-452°F). SPECIAL CARE MUST BE TAKEN DURING THE TRANSFILLING OF THE MAGNET WITH CRYOGENS AND NORMAL EXHAUST OF CRYOGENS FROM THE SYSTEM. ASIDE FROM THE OBVIOUS DANGER OF FREEZING, HELIUM GAS WILL ALSO DISPLACE THE OXYGEN IN THE ROOM. THE INSTALLATION OF AN APPROVED TOXGARD OXYGEN MONITORING SYSTEM IS RECOMMENDED.
- 2) A TRANSPORT ROUTE FOR DELIVERY OF CRYOGENS TO THE EXAM ROOM. SPECIAL VESSELS CALLED DEWARs ARE USED TO TRANSPORT HELIUM. A 250 LITER DEWAR WEIGHS 335 POUNDS AND HAS A 32" DIAMETER, A 500 LITER IS 540 POUNDS, AND 42" IN DIAMETER.
- 3) OUTSIDE VENTING OF THE HELIUM IS TO BE PROVIDED BY MEANS OF A VENT PIPE OF NON-MAGNETIC MATERIAL CALLED A QUENCH VENT.

QUENCH VENT NOTES

- 1) IN THE EVENT OF A QUENCH, THE THERMAL ENERGY DISSIPATED CAUSES AN EXTREMELY RAPID BOIL OFF OF THE LIQUID HELIUM. THE SYSTEM MUST BE CAPABLE OF VENTING THE LARGE VOLUME OF GAS GENERATED AT THE APPROXIMATE EXPANSION RATIO OF 1:700 FROM LIQUID AT 4.2°K TO ROOM TEMPERATURE GAS. THE EXHAUST SYSTEM IS CRITICAL FOR THE SAFE OPERATION OF THE MAGNET, THE DATA IN THIS DOCUMENT MUST BE FOLLOWED. SINCE HELIUM VENTED IN A QUENCH IS AN ASPHYXIANT & AN EXTREMELY COLD GAS, THE QUENCH TUBE MUST ALWAYS END AT A POINT WHERE ACCESS BY PEOPLE IS NOT POSSIBLE. QUENCH TUBE PLANNING MUST ONLY BE DONE BY QUALIFIED PERSONNEL. IT IS THE OWNER'S RESPONSIBILITY TO ENSURE THAT THE QUENCH TUBE IS MAINTAINED IN AN OPERABLE STATE.
- 2) IF THE QUENCH VENT IS NOT CONFIGURED CORRECTLY THERE IS A RISK OF DANGER THAT MAY LEAD TO DEATH OR SERIOUS INJURY AND CAN RESULT IN STRUCTURAL DAMAGE. THE EXHAUST MUST NOT BE VENTED IN AN ENCLOSED SPACE. THE OPERATOR OF THE SYSTEM MUST PREPARE AN EMERGENCY PLAN IN THE EVENT OF A QUENCH.
- 3) THE QUENCH TUBE CONSISTS OF STRAIGHT, HYDRAULICALLY SMOOTH SECTIONS, BENDS UP TO 90° AND A DIFFUSER, IF REQUIRED. THE END OF THE TUBE MUST BE PROTECTED FROM RAIN, SNOW, AND FOREIGN OBJECTS.
- 4) THE SIEMENS MAGNET HAS A QUENCH VALVE ASSEMBLY FOR CONNECTION TO THE TUBE LOCATED AT THE TOP LEFT SIDE OF THE MAGNET (SEE MAGNET ELEVATION). THE MECHANICAL CONTRACTOR WILL SUPPLY AND INSTALL A QUENCH VENT TUBE WITH CAP, TO BE NON-MAGNETIC STAINLESS STEEL (>22 GAUGE RECOMMENDED), GRADES A13304, 309, 316, OR 321 ONLY. THERMAL CONDITIONS MAY CAUSE THE TUBE TO CONTRACT UP TO 3mm/METER SO A STAINLESS STEEL BELLOWS OR FLEXIBLE SECTION MUST BE INSTALLED A MINIMUM OF EVERY 32'-9" NOT TO EXCEED 2% OF THE OVERALL LENGTH. THE MOVEMENT OF THE BELLOWS MUST BE RESTRICTED TO PREVENT EXCESSIVE EXPANSION DUE TO PRESSURE. THE WEIGHT OF THE TUBE MUST BE SUPPORTED BY THE BUILDING AND BE FLEXIBLE ENOUGH TO ALLOW MOVEMENT FROM THERMAL CONTRACTION. THE WALL EXIT SHOULD ALSO BE FLEXIBLE.
- 5) THE MAXIMUM INTERNAL PRESSURE IS CALCULATED AT 1.45 PSI. THE MAXIMUM PRESSURE SHOULD BE ENGINEERED FOR 6.5 PSI.
- 6) USE THE QUENCH VENT CALCULATOR PROVIDED BY SIEMENS TO DESIGN A QUENCH VENT THAT MEETS DESIGN REQUIREMENTS FOR DIAMETER, LENGTH, NUMBER OF ELBOWS AND PRESSURE DROP. ALL BENDS MUST BE SMOOTH WALLED AND HAVE A CENTERLINE TO INTERNAL PIPE DIAMETER RATIO OF 1.5 TO 5.0. EXPANSIONS TO PIPE DIAMETER CAN BE DONE WITH A DIFFUSER. ONLY ROUND TUBE SECTIONS MAY BE USED, RECTANGULAR SECTIONS ARE NOT ALLOWED.
- 7) THERE MUST BE A 12-19 INCH FLEXIBLE SECTION OF PIPE FOR CONNECTION TO THE QUENCH VALVE AT THE MAGNET WITH AN INSIDE DIAMETER GREATER THAN 4" AND ABLE TO WITHSTAND 6.5 PSI.
- 8) SECTIONS OF THE PIPE CAN ONLY BE JOINED BY WELDING OR BOLTED FLANGES WITH FIBER GASKETS. ROTARY FLANGES ARE PERMITTED, VEE CLAMPED FLANGES MAY NOT BE USED.
- 9) THE PROTECTION AT THE END OF THE TUBE SHALL BE 3/8" WIRE MESH COVERING AN AREA AT LEAST 2.5 TIMES THE CROSS SECTION AREA OF THE QUENCH PIPE.
- 10) WHERE THE QUENCH TUBE EXITS THROUGH A FLAT ROOF, THE THE OUTLET MUST BE ABOVE A LEVEL WHERE WATER COULD ENTER IN THE EVENT THAT THE ROOF DRAINS BECOME BLOCKED. IN THE CASE OF A HORIZONTAL EXIT THROUGH A WALL, THE OUTLET SHOULD BE ANGLED DOWNWARD NOT LESS THAN 1 PIPE DIAMETER TO PREVENT RAIN INGRESS. THE EXIT SHOULD BE LOCATED ABOVE THE LEVEL OF DRIFTING SNOW.
- 11) WHERE THE QUENCH TUBE EXITS VERTICALLY, A RAIN COVER MUST ALSO BE FITTED WITH THE DIAMETER TO BE TWO TIMES THE DIAMETER OF THE QUENCH TUBE. THE CLEARANCE BETWEEN THE RAIN GUARD AND THE MESH SHOULD 2 TIMES THE DIAMETER OF THE TUBE. A DEFLECTOR PLATE SHOULD BE WELDED TO THE TUBE WHERE IT EXITS THE ROOF TO PREVENT HELIUM FROM RE-ENTERING THE BUILDING. THE DEFLECTOR SHOULD BE AT LEAST THE DIAMETER OF THE RAIN GAURD AND LOCATED TWO PIPE DIAMETERS ABOVE THE ROOF AND TWO PIPE DIAMETERS BELOW THE RAIN GUARD.
- 12) TO AVOID INJURY FROM COLD BURNS AND ASPHYXIATION ACCESS TO THE QUENCH VENT MUST BE RESTRICTED BY 9'-11" ON EACH SIDE AND BELOW, AND 19'-9" ABOVE WITH WARNING SIGNS. THE EXIT MUST NOT LOCATED WHERE HELIUM GAS COULD BE DRAWN INTO AN AIR INLET OR OPEN WINDOW.
- 13) THE QUENCH TUBE MUST HAVE MINIMUM 1" INSULATION FOR THE FULL LENGTH. WITHIN THE RF ROOM THERE SHOULD BE A 1" LAYER OF MINERAL FIBER INSULATION WITH A VAPOR BARRIER AND 1" CLASS O OR CLASS AP ARMAFLEX. OUTDOOR PIPES MUST BE WEATHERPROOF. THE TUBE MUST HAVE A WARNING POSTED ALONG IT'S ENTIRE LENGTH FOR EXTREMELY COLD HELIUM GAS - AUTHORIZED PERSONNEL ONLY. THE INSULATION MUST NOT TOUCH THE MAGNET COVERS.
- 14) GALVANIC SEPARATION MUST BE PROVIDED BETWEEN THE MR SYSTEM, THE RF ROOM, AND THE BUILDING, TWO SEPARATIONS ARE REQUIRED USING STAINLESS STEEL BOLTS, INSULATING BUSHES AND LOCKING NUTS. NO OTHER DESIGNS ARE PERMITTED FOR SAFETY.
- 15) THE DESIGN AND CONSTRUCTION OF THE QUENCH PIPE MUST BE DOCUMENTED WITH DRAWINGS AND CALCULATIONS THAT ARE KEPT WITH INSTALLATION DOCUMENTS. IT MUST COMPLY WITH THE REQUIREMENTS IN THIS DOCUMENT BEFORE BEING CONNECTED TO THE MAGNET.

THE QUENCH VENT CONSISTS OF STRAIGHT SECTIONS AND ELBOW OF A DETERMINED DIAMETER AND LENGTH, FROM THE MAGNET QUENCH VALVE TO THE EXHAUST OUTSIDE THE BUILDING. IT SHALL BE DESIGNED BY THE MECHANICAL CONTRACTOR TO WITHSTAND A MAXIMUM PRESSURE OF 6.5 PSI (SEE QUENCH VENT NOTES-THIS SHEET). THE QUENCH VENT DESIGN MUST PASS THE QUENCH VENT CALCULATOR, PROVIDED BY THE SIEMENS PROJECT MANAGER. A SUCCESSFUL DESIGN WILL COMBINE THE LENGTH OF TUBE, DIAMETER OF TUBE, AND NUMBER OF ELBOWS SO THAT THE MAXIMUM PRESSURE IS NOT EXCEEDED.



QUENCH VENT DESIGN

NO SCALE


THE EQUIPMENT CONFIGURATION SHOWN ON THESE ELEVATIONS REPRESENT "TYPICAL" INSTALLATION CONDITIONS.

THE QUENCH VENT IS TO BE SUPPLIED AND INSTALLED BY THE MECHANICAL CONTRACTOR. IT MAY BE CONNECTED TO THE FLANGE AT THE MAGNET HORIZONTALLY OR VERTICALLY AND SHOULD BE DESIGNED, CONSTRUCTED AND INSTALLED AS PER THE INSTRUCTIONS ON THIS SHEET.

FLEXIBLE BELLOWS/ACOUSTIC DECOUPLER SECTION, FLEX COUPLING AT MAGNET PROVIDED BY SIEMENS. 90 DEGREE VERTICAL ELBOW IS PART OF STANDARD DELIVERY, HORIZONTAL ADAPTER REQUIRES SPECIAL ORDER #10140255.

14 5/8" 24 5/16" HORIZONTAL QUENCH VENT CONNECTION, 20 3/8" LEFT OF ISOCENTER

16 1/4" VERTICAL QUENCH VENT CONNECTION, 25 9/16" LEFT OF ISOCENTER



THE LOCATION, PITCH, AND MOUNTING HEIGHT ABOVE FINISHED FLOOR FOR THE MECHANICAL SYSTEMS SHALL BE SPECIFIED, DETAILED AND NOTED BY THE MECHANICAL ENGINEER OF RECORD. ALL MECHANICAL SYSTEM LOCATIONS SHALL BE COORDINATED WITH THE LOCATION OF THE CABLE TRAYS AS SHOWN ON THE 1/4" SCALE ELECTRICAL PLAN.


NOTE: COMPRESSOR LINES RUNNING FROM THE COMPRESSOR OR SEP CABINET TO THE MAGNET COLD HEAD SHALL BE ELECTRICALLY ISOLATED WITH 1/2" ARMAFLEX OR EQUAL, SUPPLIED BY CONTRACTOR AND INSTALLED BY SIEMENS. LENGTH REQUIRED SHALL BE COORDINATED WITH SIEMENS FIELD SERVICE ENGINEER AT TIME OF EQUIPMENT INSTALLATION.

MAGNET SIDE ELEVATION

SCALE: NONE

HELIUM CONTENT

MAGNET TYPE	OR 122	
LITERS AT 100%	980	=270 POUNDS
TYPICAL BOIL OFF RATE	0.01 L/HR	FOR TYPICAL CLINICAL USE, DEPENDING ON SEQUENCES AND OPERATING TIME.
TYPICAL REFILL INTERVAL	10 YEARS	




IN THE CASE OF A HORIZONTAL EXIT THROUGH A WALL, THE OUTLET MUST BE TURNED DOWN BY NOT LESS THAN THE LINE OF THE DIAMETER TO PREVENT RAIN INGRESS. THE EXIT MUST BE SITUATED WHERE IT CANNOT BE BLOCKED BY DRIFTING SNOW. TO AVOID RISK OF INJURY FROM COLD BURNS AND ASPHYXIATION, ACCESS TO THE QUENCH VENT OUTLET MUST BE RESTRICTED AS SHOWN. THE OUTLET MUST NOT BE SITUATED WHERE, IN THE EVENT OF A QUENCH, HELIUM GAS COULD BE DRAWN INTO AN AIR INLET OR OPEN WINDOW. WHERE WINDOWS ARE WITHIN THE RESTRICTED AREA, THEY MUST BE PERMANENTLY CLOSED.

OUTLET SAFETY CLEARANCES


QUENCH VENT

NO SCALE



HORIZONTAL OUTLET

SCALE: NONE



EXPANSION TO A GREATER DIAMETER TUBE CAN BE ACCOMPLISHED BY USE OF A DIFFUSER. THE DIFFUSER MUST MEET QUENCH VENT REQUIREMENTS AND THE GEOMETRIC PARAMETERS SHOWN HERE. THE DIAMETER OF THE PIPE SHOULD NEVER DECREASE IN THE DIRECTION OF FLOW.

QUENCH VENT DIFFUSER

NO SCALE

CUSTOMER/CONTRACTOR FURNISHED AND INSTALLED CRYOGEN EXITING DUCT SYSTEM.

NOTE: THE SIZE OF THE CRYOGEN EXITING DUCT MAY VARY IN SIZE DEPENDING ON THE LENGTH OF RUN AND THE AMOUNT OF BENDS REQUIRED TO EXIT THE BUILDING. SEE THE QUENCH VENT ROUTING DETAIL THIS SHEET FOR INFORMATION ON SIZING THE CRYOGEN EXITING DUCT.

QUENCH VENT FLEXIBLE SECTION TO BE DESIGNED IN ACCORDANCE WITH QUENCH VENT NOTES, THIS SHEET, BUILT BY THE MECHANICAL CONTRACTOR AND INSTALLED BY THE CUSTOMER/ CONTRACTOR. FINAL CONNECTION OF THE CRYOGEN EXITING DUCT SYSTEM SHALL BE PERFORMED BY THE CUSTOMER/CONTRACTOR DURING THE DAY OF MAGNET DELIVERY.

DIELECTRIC ISOLATION SEE DETAIL A-502/2.

RF WAVEGUIDE MADE WITH FLANGES FOR QUENCH VENT COUPLING TO BE SIZED, FABRICATED AND INSTALLED BY RF SHIELD CONTRACTOR IN COOPERATION WITH QUENCH VENT MANUFACTURER.

RF SHIELD

FLANGE FOR QUENCH VENT COUPLING

FLEXIBLE BELLOWS/ACOUSTIC DECOUPLER SECTION, FLEX COUPLING AT MAGNET PROVIDED BY SIEMENS.

MOUNTING HARDWARE TO BE SUPPLIED BY SIEMENS AND INSTALLED BY THE CUSTOMER/ CONTRACTOR. TYPICAL OF 12 PLACES. GALVANIC ISOLATION AT THIS LOCATION.

ELEVATION - SEE DETAIL 3


SIEMENS SUPPLIED AND INSTALLED ELBOW.

CRYOGEN EXITING DUCT SYSTEM NOTES:

- 1) THE MECHANICAL ENGINEER OF RECORD SHALL BE RESPONSIBLE FOR THE DESIGN, DETAIL AND SPECIFICATION OF THE CRYOGEN EXITING DUCT SYSTEM WITH ALL WORK TO PERFORMED BY THE CUSTOMER/CONTRACTOR UNLESS SPECIFIED OTHERWISE.
- 2) TAKE PROPER PRECAUTIONS WHEN FIELD WELDING IN THE VICINITY OF LIQUIFIED/FROZEN OXYGEN.
- 3) UNLESS SPECIFIED OTHERWISE ALL MATERIAL AND WORK SHALL BE PERFORMED BY THE CUSTOMER/CONTRACTOR WITH FINAL CONNECTION OF THE CRYOGEN EXITING DUCT SYSTEM TO THE SIEMENS MAGNET TO BE PERFORMED BY THE CUSTOMER/CONTRACTOR UNDER SIEMENS' SUPERVISION.
- 4) THE DIAMETER OF THE CRYOGEN EXITING DUCT MAY VARY WITH ACTUAL ROUTING.

CRYOGEN EXITING DUCT SYSTEM DETAIL

SCALE: NO SCALE



• CONNECTION OF THE QUENCH VENT PIPE TO FLANGE TO BE DESIGNED, DETAILED, AND SPECIFIED BY THE MECHANICAL ENGINEER OF RECORD TO BE FABRICATED AND INSTALLED BY THE MECHANICAL CONTRACTOR UNDER SIEMENS SUPERVISION. THE 90° ELBOW IS PART OF THE DELIVERY VOLUME, THE TWO FLANGES ARE IDENTICAL.

CONNECTING FLANGE

QUENCH VENT

NO SCALE

ATTENTION:

— THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

— THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

— IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

— ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.

— THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

SYM	DATE	DESCRIPTION
△		
—ISSUE BLOCK—		

SIEMENS
MAGNETOM ESPREE
TYPICAL FINAL DRAWINGS
MRI SUITE

THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS' AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.
ALL RIGHTS ARE RESERVED.

PROJECT #:
04103

SHEET:
M-501

DATE: 03/17/11
CHECKED: B. HERRMANN

SCALE: AS NOTED
REF. #: ---

ESPREE
03/01/11